

NATIONAL STANDARD METHOD

ENUMERATION OF ENTEROCOCCI BY MEMBRANE FILTRATION

W 3

Issued by Standards Unit, Evaluations and Standards Laboratory
Centre for Infections

ENUMERATION OF ENTEROCOCCI BY MEMBRANE FILTRATION

Issue no: 3 Issue date: 11.08.06 Issued by Standards Unit, Evaluations and Standards Laboratory in conjunction with the
Regional Food, Water and Environmental Coordinators Forum Page 1 of 13

Reference no: W 3i3

This SOP should be used in conjunction with the series of SOPs from the Health Protection Agency

www.evaluations-standards.org.uk

Email: standards@hpa.org.uk

STATUS OF NATIONAL STANDARD METHODS

National Standard Methods, which include standard operating procedures (SOPs), algorithms and guidance notes, promote high quality practices and help to assure the comparability of diagnostic information obtained in different laboratories. This in turn facilitates standardisation of surveillance underpinned by research, development and audit and promotes public health and patient confidence in their healthcare services. The methods are well referenced and represent a good minimum standard for clinical and public health microbiology. However, in using National Standard Methods, laboratories should take account of local requirements and may need to undertake additional investigations. The methods also provide a reference point for method development.

National Standard Methods are developed, reviewed and updated through an open and wide consultation process where the views of all participants are considered and the resulting documents reflect the majority agreement of contributors.

Representatives of several professional organisations, including those whose logos appear on the front cover, are members of the working groups which develop National Standard Methods. Inclusion of an organisation's logo on the front cover implies support for the objectives and process of preparing standard methods. The representatives participate in the development of the National Standard Methods but their views are not necessarily those of the entire organisation of which they are a member. The current list of participating organisations can be obtained by emailing standards@hpa.org.uk.

The performance of standard methods depends on the quality of reagents, equipment, commercial and in-house test procedures. Laboratories should ensure that these have been validated and shown to be fit for purpose. Internal and external quality assurance procedures should also be in place.

Whereas every care has been taken in the preparation of this publication, the Health Protection Agency or any supporting organisation cannot be responsible for the accuracy of any statement or representation made or the consequences arising from the use of or alteration to any information contained in it. These procedures are intended solely as a general resource for practising professionals in the field, operating in the UK, and specialist advice should be obtained where necessary. If you make any changes to this publication, it must be made clear where changes have been made to the original document. The Health Protection Agency (HPA) should at all times be acknowledged.

The HPA is an independent organisation dedicated to protecting people's health. It brings together the expertise formerly in a number of official organisations. More information about the HPA can be found at www.hpa.org.uk.

The HPA aims to be a fully Caldicott compliant organisation. It seeks to take every possible precaution to prevent unauthorised disclosure of patient details and to ensure that patient-related records are kept under secure conditions¹.

More details can be found on the website at www.evaluations-standards.org.uk. Contributions to the development of the documents can be made by contacting standards@hpa.org.uk.

Please note the references are now formatted using Reference Manager software. If you alter or delete text without Reference Manager installed on your computer, the references will not be updated automatically.

Suggested citation for this document:

Health Protection Agency (2006). *Enumeration of Enterococci by Membrane Filtration*. National Standard Method W 3 Issue 3. http://www.hpa-standardmethods.org.uk/pdf_sops.asp.

ENUMERATION OF ENTEROCOCCI BY MEMBRANE FILTRATION

Issue no: 3 Issue date: 11.08.06 Issued by Standards Unit, Evaluations and Standards Laboratory in conjunction with the Regional Food, Water and Environmental Coordinators Forum Page 2 of 13

Reference no: W 3i3

This SOP should be used in conjunction with the series of SOPs from the Health Protection Agency

www.evaluations-standards.org.uk

Email: standards@hpa.org.uk

INDEX

STATUS OF NATIONAL STANDARD METHODS	2
INDEX.....	3
AMENDMENT PROCEDURE	4
SCOPE OF DOCUMENT	5
INTRODUCTION	5
BACKGROUND	5
1 DEFINITIONS.....	6
2 PRINCIPLE	6
3 SAFETY CONSIDERATIONS	6
3.1 SAMPLE COLLECTION.....	6
3.2 SAMPLE TRANSPORT AND STORAGE	6
3.3 SAMPLE PROCESSING	6
4 EQUIPMENT	7
5 CULTURE MEDIA AND REAGENTS.....	7
6 SAMPLE PROCESSING	8
6.1 SAMPLE PREPARATION	8
6.2 FILTRATION AND INCUBATION	8
6.3 COUNTING OF COLONIES.....	8
6.4 CONFIRMATORY TESTS	8
7 CALCULATION OF RESULTS	9
8 REPORTING	9
9 QUALITY CONTROL.....	9
10 REFERENCE FACILITIES.....	10
11 ACKNOWLEDGEMENTS AND CONTACTS	10
APPENDIX 1: FLOWCHART SHOWING THE ENUMERATION OF ENTEROCOCCI BY MEMBRANE FILTRATION	11
REFERENCES	12

ENUMERATION OF ENTEROCOCCI BY MEMBRANE FILTRATION

Issue no: 3 Issue date: 11.08.06 Issued by Standards Unit, Evaluations and Standards Laboratory in conjunction with the
Regional Food, Water and Environmental Coordinators Forum Page 3 of 13

Reference no: W 3i3

This SOP should be used in conjunction with the series of SOPs from the Health Protection Agency

www.evaluations-standards.org.uk

Email: standards@hpa.org.uk

AMENDMENT PROCEDURE

Controlled document reference	W 3
Controlled document title	Standard Operating Procedure for Enumeration of Enterococci by Membrane Filtration

Each National Standard Method has an individual record of amendments. The current amendments are listed on this page. The amendment history is available from standards@hpa.org.uk.

On issue of revised or new pages each controlled document should be updated by the copyholder in the laboratory.

Amendment Number/ Date	Issue no. Discarded	Insert Issue no.	Page	Section(s) involved	Amendment
5/ 11.08.06	2.2	3	5	Background	Updated
			6	Definitions	Updated
			8 & 9	6.4.1 Aesculin hydrolysis	Updated in line with SCA 2006 Photographs of aesculin hydrolysis on Bile aesculin agar added
			12	References	Updated and incorporated into Reference Manager Software

ENUMERATION OF ENTEROCOCCI BY MEMBRANE FILTRATION

Issue no: 3 Issue date: 11.08.06 Issued by Standards Unit, Evaluations and Standards Laboratory in conjunction with the Regional Food, Water and Environmental Coordinators Forum Page 4 of 13

Reference no: W 3i3

This SOP should be used in conjunction with the series of SOPs from the Health Protection Agency

www.evaluations-standards.org.uk

Email: standards@hpa.org.uk

STANDARD OPERATING PROCEDURE FOR THE ENUMERATION OF ENTEROCOCCI BY MEMBRANE FILTRATION

SCOPE OF DOCUMENT

The method gives general guidelines for the enumeration of enterococci in samples from water intended for human consumption by the membrane filtration technique. It is also applicable to environmental and bathing waters. See QSOP 57² for examples of samples which can be examined.

INTRODUCTION

Background

Enterococci do not multiply in water. They are usually more resistant than *Escherichia coli* to environmental stress and chlorination, and as a result survive longer. Enterococci normally inhabit the intestine of man and animals, and can provide an indication of past pollution. Examination for enterococci also assists the interpretation of doubtful results from other tests such as the occurrence of large numbers of coliforms in the absence of *E. coli*³.

For potable waters, enterococci should not be detected in 250 mL samples of bottled water⁴, water in containers⁵ and natural mineral water⁶ or 100 mL of other potable waters⁷⁻⁹. Guidelines are given for bathing waters in The Bathing Waters (Classification) Regulations 1991¹⁰ and the revised Bathing Water Directive 2006/07/EC¹¹. This method is based on the method described in The Microbiology of Drinking Water, 2006³ and BS EN ISO 7899-2¹² which is the method specified in the EC Directive¹³ and the Water Supply Regulations¹⁴.

The Microbiology of Drinking Water, 2006³ contains validation data from an inter-laboratory trial for the membrane filter transfer technique for the confirmation of enterococci isolated from various waters. Laboratories should perform an in-house verification demonstrating that the method is fit for purpose.

ENUMERATION OF ENTEROCOCCI BY MEMBRANE FILTRATION

Issue no: 3 Issue date: 11.08.06 Issued by Standards Unit, Evaluations and Standards Laboratory in conjunction with the
Regional Food, Water and Environmental Coordinators Forum Page 5 of 13

Reference no: W 3i3

This SOP should be used in conjunction with the series of SOPs from the Health Protection Agency

www.evaluations-standards.org.uk

Email: standards@hpa.org.uk

1 DEFINITIONS

For the purpose of this method the following definitions apply:

Presumptive enterococci

Organisms capable of aerobic growth and of forming characteristic pink to maroon coloured colonies by reduction of tetrazolium chloride to the red dye formazan on a selective isolation medium containing sodium azide and glucose, at 37°C or 44°C within 48 hours. Some strains may produce pale colonies.

Note: Incubating at the higher temperature is a compromise between selectivity and sensitivity. Selectivity is better at 44°C although lower counts of enterococci may be obtained as some strains of enterococci do not grow at 44°C.

Enterococci

Presumptive enterococci that hydrolyse aesculin at 44°C in the presence of bile salts.

2 PRINCIPLE

A measured volume of the sample or a dilution is filtered through a membrane that is capable of retaining enterococci. The membrane is incubated on a selective/differential agar and characteristic colonies are counted. Confirmatory tests are then carried out and the result is calculated as the colony count per 100 mL of sample. For bottled waters, water in containers and natural mineral water the count is calculated as the colony count per 250 mL of sample.

3 SAFETY CONSIDERATIONS¹⁵⁻²⁴

Normal microbiology laboratory precautions apply.

3.1 SAMPLE COLLECTION

N/A

3.2 SAMPLE TRANSPORT AND STORAGE

Compliance with current transportation regulations is essential.

3.3 SAMPLE PROCESSING

- Care must be taken when using boiling water or a boiling waterbath for disinfection of the filtration apparatus
- Use heat and water resistant protective gloves and do not put the face or hands over the bath when opening to remove objects
- In addition, membrane enterococcus agar (mEA) (Section 5) contains sodium azide. This substance is highly toxic. However this medium contains only 0.04% azide and therefore has low toxicity. Although sodium azide can form spontaneously explosive compounds when it comes into contact with metals such as in copper pipes the concentration is so little that handling and disposal of the prepared medium presents no significant risk

The above guidance should be supplemented with local COSHH and risk assessments

ENUMERATION OF ENTEROCOCCI BY MEMBRANE FILTRATION

Issue no: 3 Issue date: 11.08.06 Issued by Standards Unit, Evaluations and Standards Laboratory in conjunction with the Regional Food, Water and Environmental Coordinators Forum Page 6 of 13

Reference no: W 3i3

This SOP should be used in conjunction with the series of SOPs from the Health Protection Agency

www.evaluations-standards.org.uk

Email: standards@hpa.org.uk

4 EQUIPMENT

Usual laboratory equipment and in addition:

- Membrane filtration manifold
- Filter funnels graduated to 50 mL and 100 mL
- Pyrex vacuum flask with protective jacket or equivalent: large volume eg 5 L
- Vacuum pump with moisture trap or protective filter, or alternative vacuum source
- Stainless steel flat tipped forceps
- Boiling waterbath (for instrument disinfection)
- Cycling incubator: 37°C ± 1.0°C and 44°C ± 0.5°C (optional)
- Incubator: 37°C ± 1.0°C
- Incubator: 44°C ± 0.5°C
- Petri dishes
- Cellulose ester 0.45 µm gridded membrane filters
- Automatic pipettors and associated sterile tips capable of delivering up to 10 mL and 1 mL volumes (optional)
- Pipettes, (sterile total delivery) 10 mL and 1 mL graduated in 0.1 mL volumes (optional)

5 CULTURE MEDIA AND REAGENTS

Equivalent commercial dehydrated media may be used; follow the manufacturer's instructions.

Peptone saline diluent (Maximum recovery diluent)

Peptone	1.0 g
Sodium chloride (NaCl)	8.5 g
Water	1 L
pH	7.0 ± 0.2 at 25°C

Quarter strength Ringer's solution

Sodium chloride (NaCl)	2.25 g
Potassium chloride (KCl)	0.11 g
Calcium chloride (CaCl ₂)	0.12 g
Sodium bicarbonate (NaHCO ₃)	0.05 g
Water	1 L

Membrane enterococcus agar (mEA)

Tryptose	20.0 g
Yeast extract	5.0 g
Glucose	2.0 g
Dipotassium hydrogen phosphate (K ₂ HPO ₄)	4.0 g
Sodium azide (NaN ₃)	0.4 g
2,3,5 triphenyl tetrazolium chloride (TTC)	0.1 g
Agar	12.0 g
Water	1 L
pH	7.2 ± 0.2 at 25°C

Bile aesculin agar (BAA)

Bile salts	20.0 g
Peptone	8.0 g
Ammonium iron (III) citrate	0.5 g
Aesculin	1.0 g

ENUMERATION OF ENTEROCOCCI BY MEMBRANE FILTRATION

Issue no: 3 Issue date: 11.08.06 Issued by Standards Unit, Evaluations and Standards Laboratory in conjunction with the Regional Food, Water and Environmental Coordinators Forum Page 7 of 13

Reference no: W 3i3

This SOP should be used in conjunction with the series of SOPs from the Health Protection Agency

www.evaluations-standards.org.uk

Email: standards@hpa.org.uk

Agar	12.0 g
Water	1 L
pH	7.1 ± 0.2 at 25°C

6 SAMPLE PROCESSING

6.1 SAMPLE PREPARATION

Water samples are to be received and handled as described in SOP W 1 Section 5. In brief the nature of the request and condition of the sample is to be noted on arrival. Protect the samples from direct sunlight and transport in an insulated container or refrigerator at 2°C – 8°C. Samples should be analysed as soon as is practicable on the day of collection. In exceptional circumstances, if there is a delay, storage under the above conditions should not exceed 24 hours before the commencement of analysis.

Following the procedures laid down in SOP W 1 Section 5 select suitable volumes for analysis and prepare any necessary dilutions.

6.2 FILTRATION AND INCUBATION

Following the procedures laid down in SOP W 1 Section 5, filter a measured volume of sample through the membrane. For bottled waters, water in containers and natural mineral waters filter 250 mL water, for other potable quality waters use volumes of 100 mL. For surface waters use 100 mL and, if expected to be highly polluted, 10 mL and 1 mL volumes.

Place the membrane onto mEA.

At 37°C ± 1°C for 44 hours ± 4 hours for potable water, if few enterococci are expected such as in drinking water.

At 37°C ± 1°C for 4 hours ± 0.5 hour followed by 40 ± 4 hours at 44°C ± 0.5°C for non-potable water likely to contain a high background count eg river water, sea water.

6.3 COUNTING OF COLONIES

After incubation, enumerate the presumptive enterococci by counting all red, maroon, pink or pale pink colonies.

6.4 CONFIRMATORY TESTS

6.4.1 AESCULIN HYDROLYSIS

Transfer the membrane filter from mEA (membrane filter face up with the colonies uppermost) to a BAA plate pre-warmed to room temperature. Ensure that no air bubbles are trapped between the membrane filter and the medium. Incubate the BAA plate at 44°C ± 0.5°C for 4 hrs ± 5 minutes.

ENUMERATION OF ENTEROCOCCI BY MEMBRANE FILTRATION

Issue no: 3 Issue date: 11.08.06 Issued by Standards Unit, Evaluations and Standards Laboratory in conjunction with the Regional Food, Water and Environmental Coordinators Forum Page 8 of 13


Reference no: W 3i3

This SOP should be used in conjunction with the series of SOPs from the Health Protection Agency

www.evaluations-standards.org.uk

Email: standards@hpa.org.uk

Figure 1 Colonies on bile aesculin agar (BAA)



Enterococcus faecalis (large colonies with halo) with *Aerococcus viridans* (small colonies without halo)

Development of colonies that are surrounded by a black or brown halo in the surrounding medium resulting from the aesculin hydrolysis process confirms the presence of enterococci (see Figure 1). The use of a heavy inoculum will increase the speed at which results are obtained, and a positive result may be evident after a few hours incubation.

7 CALCULATION OF RESULTS

Calculate the presumptive count of enterococci as follows:

$$\text{Presumptive count / 100mL} = \frac{\text{Number of colonies counted}}{\text{Volume tested}} \times 100$$

Following the procedures described in SOP W 1 Section 7, calculate the number of confirmed enterococci present in the original sample.

8 REPORTING

Report the results using the procedure described in SOP W 1 Section 9.

If enterococci are not detected, report as:

Not detected in 100 mL

If the test organisms are detected report as:

a in 100 mL

where **a** is the presumptive or confirmed count

For natural mineral waters the result is reported as the count per 250 mL.

9 QUALITY CONTROL

Membrane filtration

When the membrane filtration technique is used internal quality control procedures must be carried out at least once daily depending on the workload of the laboratory. If more than one batch of media is used in a session it is necessary to repeat the quality control test for each batch.

ENUMERATION OF ENTEROCOCCI BY MEMBRANE FILTRATION

Issue no: 3 Issue date: 11.08.06 Issued by Standards Unit, Evaluations and Standards Laboratory in conjunction with the Regional Food, Water and Environmental Coordinators Forum Page 9 of 13

Reference no: W 3i3

This SOP should be used in conjunction with the series of SOPs from the Health Protection Agency

www.evaluations-standards.org.uk

Email: standards@hpa.org.uk

The quantitative internal quality controls are to be carried out using suspensions of positive and negative control organisms known to contain less than 100 colony forming units in the volume filtered.

Positive control:

Enterococcus faecalis NCTC 775

Negative control:

Escherichia coli NCTC 9001

Blank control:

Filter 100mL of sterile distilled water, peptone saline diluent or quarter strength Ringer's solution using the same funnel as used for the positive control following disinfection.

Incubate all tests in parallel with routine tests, proceed with confirmatory tests and determine the count.

Confirmatory tests

Details of the test strains are presented below. Inoculate each batch of medium/reagent used on the day of testing with the control organism. Incubate in parallel with the routine tests.

Bile aesculin agar

Test strains 44°C ± 0.5°C for the same incubation time as the routine tests.

Enterococcus faecalis NCTC 775

positive (black)

Escherichia coli NCTC 9001

negative (colourless)

10 REFERENCE FACILITIES

N/A

11 ACKNOWLEDGEMENTS AND CONTACTS

This National Standard Method has been developed, reviewed and revised by the Water Working Group for Standard Methods (http://www.hpa-standardmethods.org.uk/wg_water.asp). The contributions of many individuals in clinical bacteriology laboratories and specialist organisations who have provided information and comment during the development of this document, and final editing by the Medical Editor are acknowledged.

The National Standard Methods are issued by Standards Unit, Evaluations and Standards Laboratory, Centre for Infections, Health Protection Agency London.

For further information please contact us at:

Standards Unit
Evaluations and Standards Laboratory
Centre for Infections
Health Protection Agency
Colindale
London
NW9 5EQ
E-mail: standards@hpa.org.uk

ENUMERATION OF ENTEROCOCCI BY MEMBRANE FILTRATION

Issue no: 3 Issue date: 11.08.06 Issued by Standards Unit, Evaluations and Standards Laboratory in conjunction with the Regional Food, Water and Environmental Coordinators Forum Page 10 of 13


Reference no: W 3i3

This SOP should be used in conjunction with the series of SOPs from the Health Protection Agency

www.evaluations-standards.org.uk

Email: standards@hpa.org.uk

APPENDIX 1: FLOWCHART SHOWING THE ENUMERATION OF ENTEROCOCCI BY MEMBRANE FILTRATION



ENUMERATION OF ENTEROCOCCI BY MEMBRANE FILTRATION

Issue no: 3 Issue date: 11.08.06 Issued by Standards Unit, Evaluations and Standards Laboratory in conjunction with the Regional Food, Water and Environmental Coordinators Forum Page 11 of 13

Reference no: W 3i3

This SOP should be used in conjunction with the series of SOPs from the Health Protection Agency

www.evaluations-standards.org.uk

Email: standards@hpa.org.uk

REFERENCES

1. Department of Health NHS Executive: The Caldicott Committee. Report on the review of patient-identifiable information. London. December 1997.
2. National Standard Method: The microbiological examination of water samples QSOP57. London: Health Protection Agency; 2005.
3. Standing Committee of Analysts. Environment Agency. The Microbiology of Drinking Water (2006). Part 5: The isolation and enumeration of enterococci by membrane filtration. London. http://www.environment-agency.gov.uk/commodata/acrobat/mdw2006_part_5_1312024.pdf.
4. Council Directive 98/83/EC: Relating to the Quality of Water Intended for Human Consumption. Official Journal of the European Communities; 1998. p. L330-32-L330/53
5. Council Directive 98/83/EC: Relating to the Quality of Water Intended for Human Consumption. Official Journal of the European Communities; 1998. p. L330-32-L330/53
6. The Natural Mineral Water, Spring Water and Bottled Drinking Water Regulations (Statutory Instrument No. 666). London: HMSO; 2003.
7. Council Directive 98/83/EC: Relating to the Quality of Water Intended for Human Consumption. Official Journal of the European Communities; 1998. p. L330-32-L330/53
8. The Water Supply (Water Quality) Regulations (Statutory Instrument No. 3184). London: HMSO; 2000.
9. The Drinking Water in Containers Regulations (Statutory Instrument No. 743). London: HMSO; 1994.
10. Great Britain. The Bathing Waters (Classification) Regulations 1991: Statutory Instrument No. 1597 . HMSO. http://www.opsi.gov.uk/si/si1991/Uksi_19911597_en_1.htm.
11. 2006/07/EC: Council Directive concerning the management of bathing water quality (Repealing Directive 76/160/EEC). L64. Official Journal of the European Communities. http://europa.eu.int/eur-lex/lex/LexUriServ/site/en/oj/2006/l_064/l_06420060304en00370051.pdf. p. 37-51.
12. ISO 7899-2: 2000. Water quality - Detection and enumeration of intestinal enterococci. Part 2: Membrane filtration method. International Organization for Standardization; 2000.
13. Council Directive 98/83/EC: Relating to the Quality of Water Intended for Human Consumption. Official Journal of the European Communities; 1998. p. L330-32-L330/53
14. The Water Supply (Water Quality) Regulations (Statutory Instrument No. 3184). London: HMSO; 2000.
15. Advisory Committee on Dangerous Pathogens. 2004 Approved List of Biological Agents. <http://www.hse.gov.uk/pubns/misc208.pdf>. p. 1-17.
16. Public Health Laboratory Service Standing Advisory Committee on Laboratory Safety. Safety Precautions: Notes for Guidance. 4th ed. London: Public Health Laboratory Service (PHLS); 1993.
17. Control of Substances Hazardous to Health Regulations 2002. General COSHH. Approved Code of Practice and Guidance, L5. Suffolk: HSE Books; 2002.
18. Health and Safety Executive. 5 steps to risk assessment: a step by step guide to a safer and healthier workplace, IND (G) 163 (REVL). Suffolk: HSE Books; 2002.

ENUMERATION OF ENTEROCOCCI BY MEMBRANE FILTRATION

Issue no: 3 Issue date: 11.08.06 Issued by Standards Unit, Evaluations and Standards Laboratory in conjunction with the Regional Food, Water and Environmental Coordinators Forum Page 12 of 13

Reference no: W 3i3

This SOP should be used in conjunction with the series of SOPs from the Health Protection Agency

www.evaluations-standards.org.uk

Email: standards@hpa.org.uk

19. Health and Safety Executive. A guide to risk assessment requirements: common provisions in health and safety law, IND (G) 218 (L). Suffolk: HSE Books; 2002.
20. NHS Estates. Health Building Note 15. Facilities for pathology services. 2nd ed. London: Her Majesty's Stationary Office (HMSO); 2005.
21. BS EN 12469: 2000. Biotechnology - performance criteria for microbiological safety cabinets. London: British Standards Institution (BSI); 2000.
22. BS 5726: 1992. Microbiological safety cabinets. Part 2. Recommendations for information to be exchanged between purchaser, vendor and installer and recommendations for installation. London: British Standards Institution (BSI); 1992.
23. BS 5726: 1992. Microbiological safety cabinets. Part 4. Recommendations for selection, use and maintenance. London: British Standards Institution (BSI); 1992.
24. Advisory Committee on Dangerous Pathogens. The management, design and operation of microbiological containment laboratories. Suffolk: HSE Books; 2001.

ENUMERATION OF ENTEROCOCCI BY MEMBRANE FILTRATION

Issue no: 3 Issue date: 11.08.06 Issued by Standards Unit, Evaluations and Standards Laboratory in conjunction with the Regional Food, Water and Environmental Coordinators Forum Page 13 of 13

Reference no: W 3i3

This SOP should be used in conjunction with the series of SOPs from the Health Protection Agency

www.evaluations-standards.org.uk

Email: standards@hpa.org.uk